

# Integrating Digital Solutions into Implant Therapy

## From Scan to Insertion





# Objectives

- 1 Integration to a digital workflow
- 2 Hands-on experience
- 3 Surgical & restorative instrumentation
- 4 Digital workflow outsourcing







# Stations 1 & 5

*Pre-Treatment Scanning  
& Implant Scanning*

# Station 2

*Planning*

# Station 6

*Finalization  
& Insertion*

# Station 4

*Instrumentation*

# Station 3

*Design & 3D  
Printing/Milling*



# Station I

## *Pre-Treatment Scanning*



Scanning tips



Hands on trial



Task delegation



# Station 1 Tips

## *Pre-Treatment Scanning*

### Intraoral scanning


- Sit the patient upright when scanning the bite for a more accurate occlusion
- Let the patient rinse with cool water first for less fogging of the scanner lens
- Delegate scanning to your assistants & hygienists

### CBCCT


- Take CBCCTs with an open bite to prevent overlap
- Try to minimize scatter and consider provisionalization prior to taking a CBCCT
- Take a scout film if necessary to prevent retaking the radiograph(s)

# Station 2

## Planning



CBCT merging



Implant planning



Design practice



# Station 2 Tips

## *Planning*

### Intraoral scanning

- Make sure your scans don't have holes & missing data
- Confirm the articulation is correct
- Don't rely on quadrant scans - always scan the full mouth & bite

### CBCT

- Make sure your scans are of diagnostic quality
- Consider reviewing your scans with the lab as well
- The doctors must make all final design decisions



# Station 3

## *Design & 3D Printing/Milling*



Live 3D printing



Lab perspective



Equipment demo



# Station 3 Tips

## *Design & 3D Printing/Milling*

- When planning a case, print a cast for reference to minimize chair time
- Understand what your lab needs from your intraoral scans
- Rely on your technology manufacturer's tech support
- Invest in a small 3D printer for your office if you don't have a lab onsite



# Station 4

## *Instrumentation*



Surgical Guides



Genesis Kits



Organizational Tips



# Station 4 Tips

## *Instrumentation*

### Instrumentation

- Planning for guided surgery is important - make sure your data is accurate
- Keyless system = easier guided surgery with a more predictable outcome
- Order from your stock: keep your inventory up to date

### Organization & Management

- Keep an heir and a spare of everything & everyone
- Centralize your tasks and cross-train your staff to support each other
- Get creative with your organization and label everything



# Station 5

## *Implant Scanning*



Single implant



Multiple implants



Complete arch



# Station 5 Tips

## *Implant Scanning*

- Understand scan body placement and positioning
- Consider copying existing scans to minimize chair time
- Full arch cases are more successful than ever
- Separate scan bodies for fixture level and abutment level scans
- Practice, practice, practice!




# Station 6

## *Finalization & Insertion*



Pearls for success



Materials review



Insertion demo



# Station 6 Tips

## *Finalization & Insertion*

- Ti base cementation
- Screwmentable prosthesis cementation - importance of cement removal
- Try-in and verification of implant prosthesis seating with radiographs
- Sealing of access openings with an FDA approved PTFE material
- Provisional vs permanent access closure





[gdental.com/events](http://gdental.com/events)